

# BPS Capsule Filters

PES Membrane



Validated for use in multiple biopharmaceutical applications

Excellent flow rates with high throughput

Designed for minimal leachables and extractables

Low adsorption of proteins and preservatives

## Flow Rate / Filtration Area

The following table represents typical water flow at a one psi (69 mbar) pressure differential across a single 2 inch capsule with 1.0 ft<sup>2</sup> (930 cm<sup>2</sup>) of media with 1/2" FNPT ports. The test fluid is water at ambient temperature. Higher pressure drops are acceptable, but as flows increase the pressure drop of the housing becomes more apparent.

| Pore Size | 0.03 μm | 0.10 μm | 0.22 μm | 0.45 μm | 0.65 μm | 0.8 μm | 1.0 μm | 1.2 μm |
|-----------|---------|---------|---------|---------|---------|--------|--------|--------|
| GPM       | 0.16    | 0.26    | 0.46    | 0.71    | 0.86    | 0.91   | 0.97   | 1.0    |
| LPM       | 0.61    | 0.98    | 1.74    | 2.69    | 3.26    | 3.44   | 3.67   | 3.78   |

\* For approximate flow rates for 5" through 30" capsules, refer to the appropriate cartridge data sheet

## Construction Materials

|                    |                                 |
|--------------------|---------------------------------|
| Housing            | Polypropylene                   |
| Filtration Media   | Polyethersulfone (PES) Membrane |
| Media Support      | Polypropylene                   |
| End Caps           | Polypropylene                   |
| Center Core        | Polypropylene                   |
| Outer Support Cage | Polypropylene                   |
| Sealing Method     | Thermal Bonding                 |

## Applications

- ◆ Buffers and Feedstocks
- ◆ WFI Water
- ◆ Solvents
- ◆ SVPs
- ◆ LVPs
- ◆ Vaccines

Biopharmaceutical grade PS capsules are hydrophilic and manufactured with the highest quality asymmetric polyethersulfone (PES) membrane. PES membrane exhibits excellent flow rates with high throughput. BPS capsule filters are used for bioburden reduction applications in the biologicals industry. Specific applications for BPS capsule filters are bacteria removal in buffers, USP Water for Injection (WFI), SVPs, LVPs and other biological products. The filter elements are 100% integrity tested during production.

Polyethersulfone is particularly suited for the filtration of products that contain elements that can adsorb to the media, such as preservatives and proteins. The lower binding characteristics of PES make it a good choice for filtration of valuable protein solutions such as vaccines and other biologicals.

## BPS Capsule Filters - Filtration Area

| Media        | Capsule Length                               |                                              |                                              |                                                |                                                |
|--------------|----------------------------------------------|----------------------------------------------|----------------------------------------------|------------------------------------------------|------------------------------------------------|
|              | 2"                                           | 5"                                           | 10"                                          | 20"                                            | 30"                                            |
| PES Membrane | 1.0 ft <sup>2</sup><br>(0.93m <sup>2</sup> ) | 3.0 ft <sup>2</sup><br>(2.78m <sup>2</sup> ) | 7.0 ft <sup>2</sup><br>(6.50m <sup>2</sup> ) | 14.0 ft <sup>2</sup><br>(13.01m <sup>2</sup> ) | 21.0 ft <sup>2</sup><br>(19.51m <sup>2</sup> ) |

## Sanitization/Sterilization

**Autoclave** ..... 250° F (121° C), 30 min, multiple cycles  
**Chemical Sanitization** ..... Industry standard concentrations of hydrogen peroxide, paracetic acid, sodium hypochlorite and other selected chemicals.

**Note** .....BPS capsules are not to be used in steam.  
**Pre-Sterilized** .....BPS capsules are offered in both non- and pre-sterilized forms.

## Maximum Operating Parameters

|                                |                                                |
|--------------------------------|------------------------------------------------|
| Liquid Operational Pressure    | 80 psi (5.5 bar) at 20 °C (68 °F)              |
| Gases Operational Pressure     | 60 psi (4.1 bar) at 20 °C (68 °F)              |
| Operating Temperature          | 43 °C (110 °F) at 30 psi<br>(2.1 bar) in water |
| Forward Differential Pressure  | 50 psid (3.4 bard) at 20 °C (68 °F)            |
| Reverse Differential Pressure  | 40 psid (2.7 bard) at 20 °C (68 °F)            |
| Recommended Changeout Pressure | 35 psid (2.4 bard)                             |

## Integrity Test Specifications

| Pore Size | Test Pressure (psi) | Max Diffusion Rate (cc/min -water wetted membrane) |      |     |     |     |
|-----------|---------------------|----------------------------------------------------|------|-----|-----|-----|
|           |                     | 2"                                                 | 5"   | 10" | 20" | 30" |
| 0.03      | 60                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 0.10      | 48                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 0.22      | 35                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 0.45      | 20                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 0.65      | 15                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 0.80      | 12                  | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 1.00      | 8                   | 4.3                                                | 12.9 | 30  | 60  | 90  |
| 1.20      | 7                   | 4.3                                                | 12.9 | 30  | 60  | 90  |

## Integrity Test Specifications - Bubble Point

| Pore Size | Bubble Point (water wetted membrane) |
|-----------|--------------------------------------|
| 0.03 µm   | **                                   |
| 0.10 µm   | **                                   |
| 0.22 µm   | 50 psig (3.5 barg)                   |
| 0.45 µm   | 25 psig (1.7 barg)                   |
| 0.65 µm   | 19 psig (1.3 barg)                   |
| 0.8 µm    | 15 psig (1.1 barg)                   |
| 1.0 µm    | 10 psig (0.7 barg)                   |
| 1.2 µm    | 9 psig (0.6 barg)                    |

\*\* Test pressure exceeds operational limits of capsule filters. Use the diffusion test method.

## Validation


BPS cartridges are validated using test procedures that comply with the intent of both ASTM F 838-05 and HIMA protocols for the determination of bacterial retention in filters used for liquid filtration. The challenge level is 10<sup>7</sup> organisms per cm<sup>2</sup> of filter media:

- 0.10 µm challenged with *Acholeplasma laidlawii*;
- 0.22 µm challenged with *Brevundimonas diminuta*;
- 0.45 µm challenged with *Serratia marcescens*;
- 0.65 µm challenged with *Saccharomyces cerevisiae*.

Critical Process Filtration can provide validation assistance.

## Ordering Information

Capsule number example: Biopharmaceutical Grade PES Membrane, 0.22 Micron Rating, Pre-Sterilized, 10" Length, Sanitary Inlet, Sanitary Outlet = CPBPS-20S0001FF.



Request a **QUOTE** from your area representative



### Critical Process Filtration, Inc.

One Chestnut Street • Nashua, NH 03060  
Tel: 603.880.4420 • Fax: 603.880.4536

[criticalprocess.com](http://criticalprocess.com) • [sales@criticalprocess.com](mailto:sales@criticalprocess.com)

The information contained herein is subject to change without notice.

The Critical Process Filtration logo is a trademark of Critical Process Filtration, Inc.

Viton is a trademark of DuPont Performance Elastomers L.L.C.

© 2011-2015 Critical Process Filtration, Inc. • All Rights Reserved • Data Sheet CPBPSDS1011 Rev-

## Extractables

Biopharmaceutical grade filters typically exhibit low levels of non-volatile residues.

## Quality Assurance and Standards

Critical Process Filtration uses state of the art computer controlled equipment to consistently produce high quality products as well as significantly reduce hand operations that can compromise quality. All manufacturing and testing is continuously monitored in real time so that data can be quickly and easily analyzed to facilitate improvements in both quality and cost.

The Critical Process Filtration manufacturing and quality systems meet rigorous ISO 9001:2008 standards. Each operation, including assembly, testing, cleaning, drying and packaging, is done in an appropriately rated clean room. Manufacturing is controlled using a sophisticated manufacturing system that networks work stations, manufacturing centers and inspection points. During the manufacturing and inspection processes, data is collected in real time to allow continuous quality monitoring and full traceability of all materials and processes.

Each capsule filter assembly is integrity tested before release.

## USP Biosafety and FDA Compliance

The materials used to construct biopharmaceutical grade PS capsule filters are non-toxic and meet the requirements for the MEM Elution Cytotoxicity Test and the requirements for Biological Reactivity Tests in the current version of the United States Pharmacopeia (USP) for Class VI - 121 °C Plastics. In addition, the materials meet the requirements listed by the FDA as appropriate for use in articles intended for repeated food contact as specified in Title 21 CFR sections 174.5, 177.1500, 177.1520, 177.1630, 177.2440, and 177.2600 as appropriate. BPS capsule filters comply with Title 21 CFR sections 210.3 (b)(6) and 211.72, for non-fiber releasing filters. The levels of bacterial endotoxins in aqueous extracts from biopharmaceutical grade capsule filters are below current USP limits as specified for water for injection.

## Total Performance

Critical Process Filtration, Inc. is a vertically integrated manufacturer of filtration products to industries in which filtration is considered a critical part of the manufacturing process. We supply a complete line of products and services to help you cost effectively satisfy all your filtration requirements from a single source.

### Hose Barb Diameter Ranges\*

|                        | Minimum        | Maximum         |
|------------------------|----------------|-----------------|
| <b>Outer Diameters</b> | 11/32" (8.6mm) | 9/16" (14.0mm)  |
| <b>Inner Diameters</b> | 5/32" (4.0mm)  | 13/32" (10.5mm) |